



IV Pharmacy & Nutrition Services

Charging the Power Pack

Before each use recharge the Power Pack every 24 hours for 4 hours

1. Attach the AC Adapter connector to the power pack INPUT connector. (Figure 1)
2. Plug the AC Adapter's 3-prong plug into An AC wall Outlet
3. Note the power pack indicator lights

The Amber indicator light shows the power Pack is **charging**

The green indicator light shows the power Pack is **ready** (charge for 4 hours)

Connecting the Power Pack to the CADD-Prizm Pump

1. Open the Power jack cover on the side of the pump
DO NOT USE THE DATA JACK
2. Line up the red dot on the power pack's output Connector with the red dot on the pump. Push the Connector forward until it stops. **DO NOT TWIST CONNECTOR**
3. Watch for the message "Power Pack Good" on

The pump's display

Revised 11/13/07